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Contents

1.	General instructions	12
2.	Safety instructions to be carefully followed!	12
3.	Start	
3.1.	Degassing	13
3.2.	Cleaning	13
3.3.	Ultrasound ON / OFF	15
3.4.	Units with heating	15
3.5.	Emptying the ultrasonic tank	16
4.	Fixing and removing of the outlet at TK 100 H	16
5.	Additional recommendations	
5.1.	Cleaning agents	17
5.2.	Cleaning and care	17
5.3.	Treatment of medical instruments	18
6.	Error findings ?	18
7.	Technical data	19

1. General instructions

- **Never** run the unit without liquid!
- **Do not** put parts to be cleaned into the tank, which can have aggressive soiling such as acids, chloride ions or similar substances!
- **Do not** fill the stainless steel oscillating tanks with aggressive cleaning agents such as acids or chemicals which contain or release chloride ions (some disinfectants, dish-washing detergents, household cleaners, saline solutions).
- **Do not** fill the stainless steel oscillating tanks with combustible liquids.
- Use aggressive cleaning agents **only** in inset beakers or non-perforated plastic insert tanks.
- Place **nothing** on the bottom of the oscillating tank; use accessories!
- **Do not** put your hand into the oscillating liquid!
- The ultrasonic energy in the oscillating tank heats up the cleaning liquid without additional heating! Check the temperature regularly when cleaning temperature-sensitive parts.

2. Safety instructions to be carefully followed!

- After removing the packing, check the unit's casing and control elements for any signs of damage. In such case, do not connect the unit to the electrical power supply, but report the damage immediately and in writing to the carrier and your supplier.
Save the original packing.
- Place the unit in a horizontal position on a dry, solid base.
- Plug the unit to grounded sockets only.
- Keep the unit's surface and control elements clean and dry!
- Replace defective parts only with original parts!
- Do not clean any living beings (plants and animals) in the unit!
- If preheating of the cleaning liquid is required, it is necessary to stir it at least every 30 minutes or switch on ultrasound for 5 minutes (retardation of boiling).
- Ultrasonic cleaning produces cavitation noise. If working within distances of 2 meters from the unit for longer periods of time, be sure to wear ear plugs to avoid damage to your hearing!
- Do not move or carry units filled with liquid.

Principle of ultrasonic cleaning

PZT (lead zirconate titanate) oscillators underneath the tank convert electrical energy into mechanical vibrations. units produce oscillations of 35 kHz (TK 100 /H 43 kHz) in the cleaning liquid. Tiniest vacuum bubbles are formed which then implode (cavitation). This cavitation principle removes dirt and grease from the deepest pores - "electronic scrubbing" - at all points that come into contact with the cleaning liquid.

3. Start

- Make sure that the outlet of the unit – either ball cock, closing cap (RK 100 SH) or clamp (TK 100 H) - is closed !
Fix the outlet kit at TK 100 H before starting, see point 4.
- Fill the oscillating tank up to the filling mark.
- Connect the unit to a grounded socket.

3.1. Degassing

Degassing the cleaning liquid enhances the cleaning effect.


Cleaning and contact liquids contain dissolved gases (for example, oxygen). Sonicate cleaning and contact liquid, which have just been poured into the oscillating tank or have been left in it for a longer period of time, for approx. 5 to 15 minutes before starting the cleaning process. Degas oscillating tanks with a volume of more than 10 l for roughly 30 minutes. During the degassing process, the cavitation noise changes and the loud degassing noise decreases toward the end of the process, the unit appears to run more quietly. Lower noise levels indicate the end of the degassing process, not a reduction in ultrasonic energy.

3.2. Cleaning

- Be sure that the objects in the tank are completely covered with liquid.
- Normally the objects may be cleaned **DIRECTLY** in the ultrasonic tank with an insert basket.
- In the following cases clean objects **INDIRECTLY** in inset beakers or non-perforated insert tanks to protect the stainless steel oscillating tank:
 - when chemically aggressive cleaning liquids (for example acids, etc.) are used.
 - when chemically aggressive substances (for example when cleaning racks from film developing machines) or
 - when abrasive substances (for example polishing paste, sand) are to be removed.


DIRECTLY in the oscillating tank:

- Fill the oscillating tank to 2/3 with cleaning liquid. Improperly filled units can suffer damage (see point 7 – dosing table).
- Use always an insert basket or utensil holder.
- Parts must not be placed on the tank bottom (Exemption basket K6 and SH 7).


INDIRECTLY in inset beakers or non-perforated insert plastic tub

- Fill the oscillating tank with contact liquid (for example TICKOPUR R 33 diluted with water, 1 % solution).
- Fill the inset beakers or insert plastic tub with cleaning liquid
- Suspend inset beakers in the positioning cover, do not place them on the tank bottom.
- Inset beakers or insert plastic tub must be immersed in the contact liquid to a depth of at least 2 cm. See to the tank filling to 2/3.
- Several inset beakers with different cleaning liquids can be sonicated at the same time.

When inserting the objects to be cleaned, be sure to observe the following:


- Place the dirtier side towards the bottom.
- Do not stack - ultrasonic energy will be absorbed.
- Sensitive pieces may not be allowed to touch one another. Avoid damage by using plastic insets, plastic insert baskets or the knobbed silicone mats in insert basket.
- Be sure that there are no air bubbles in hollow spaces (for example pocket holes).
- Check the temperature when cleaning heat-sensitive objects.

3.3. Ultrasound ON / OFF

Timer operation

- Turn the knob clockwise to the appropriate cleaning time → the green pilot lamp lights up.
- The timer switches off automatically.
- You can shorten the time of cleaning or switch the unit off by turning the knob anticlockwise.

Continuous operation (∞ position)

- Turn the knob anticlockwise until it clicks into position.
- The unit does not switch off automatically, turn the knob back to "0" in order to switch off.

The unit may remain plugged in when it is not running.
To disconnect it pull the mains plug.


Cleaning time

Cleaning time must be kept as short as possible. Depending on the degree of soiling and the cleaning liquid used (see the data on its label), the cleaning time varies between 1 and 15 minutes.


3.4. Units with heating

Units marked with „H“ have a built-in heating.
The heating system runs independently from the ultrasonic energy.

- Tip switch operation (fixed temperature of 65 °C):
Heating switch on – control lamp is on.
When the bath temperature of approx. 65 °C is reached, the heating will switch off, the control lamp remains on.

- Set the temperature knob clockwise to the selected temperature → the thermostatic regulation is in operation:
 - the yellow and white pilot lamps light up.
 - the yellow pilot lamp extinguishes when the selected temperature is reached.

Important:

- Warmed-up liquids improve the cleaning efficiency.
- Continuous irradiation and covering the oscillating tank can cause the temperature of the liquid rise above the set value.
- Do not heat disinfectants.
- Simultaneous irradiation shortens the heating-up period.
- The ultrasonic energy in the oscillating tank heats the cleaning liquid without any additional heating.

3.5. Emptying the ultrasonic tank


Pull the mains plug.

Do not place the unit into a sink.

Units with drain

- Ball cock

Turn the lever to the outflow direction.


- Outlet hose, closing cap or hose clamp

When unscrewing the closing cap or unlocking the hose clamp, make sure that the hose end is in a higher position than the liquid level in the oscillating tank. – It is possible, that the hose (at the clamp position) does not open immediately caused by quenching.

Units without drain


Pour out the liquid over one of the tank corners.
Avoid splashing from below.

4. Fixing and removing of the outlet at TK 100 H

TK 100 H


Back view

Outlet:

- 1 angled fitting
- 2 tube
- 3 fixing clamp
- 4 closing clamp
- 5 tube support

Detail A:

For fixing the outlet, push the angled fitting over the outlet pipe at the back of the unit. The angled fitting can be removed by pressing the small ring to the fitting with a small screw driver and pulling it from the outlet pipe, simultaneously. Do not tilt!

5. Recommendations

5.1. Cleaning agents

- Replace used liquids and disinfectants, do not "freshen them up" with an added dose.
- Layers of dirt on the tank bottom or in the insets (abrasives and polishing agents, for example, as well as other deposits) reduce the cleaning efficiency. Empty the oscillating tank and clean it regularly. Replace the cleaning liquid.
- The cleaning and contact liquids in the ultrasonic oscillating tank may not exceed a maximum operating temperature of 100 °C.
- Distilled or deionized water without additives is unsuitable for ultrasonic cleaning directly in the oscillating tank. They lack the necessary cavitation properties to form a uniform ultrasonic field. Irregular cavitation can damage the oscillating tank or the objects being cleaned. (You may use distilled or deionized water in the inset beakers and non-perforated insert tubs.)
- Use only TICKOPUR, STAMMOPUR (see above) in the oscillating tank. These agents are especially developed and recommended for ultrasonic cleaning.
- Be careful when working with aggressive cleaning agents in insets or non-perforated insert tanks; avoid splashing into the contact liquid or onto stainless steel surfaces; if necessary, replace the contact liquid immediately, clean the surfaces and wipe them dry.

The manufacturer does not take any responsibility for equipment damage caused by unsuitable cleaning chemicals.

units require no maintenance.

5.2. Cleaning and care

- Be sure to pull the plug before cleaning the unit.
- Do not shower the units or immerse them in water; instead, wipe them with a damp cloth, in- and outside.
- Rinse out the stainless steel tank frequently and wipe it dry.
- Do not use steel wool, scrapers or graters.
- If rings should appear in the stainless steel tank after a longer period of operation, remove them with a customary stainless steel cleaner without scouring additives.
- Metal parts remaining on the stainless steel surface, as well as particles of rust, can penetrate the passive protective layer of the stainless steel, thereby "activating" it and causing it to rust. This rust provokes pittings on the stainless steel surface. Therefore, remove all remaining metal parts such as screws, metal chips and the like; remove minor rust spots immediately with a soft cloth and customary stainless steel cleaners without scouring additives.

5.3. Treatment of medical instruments

- The instruments may not be placed on the tank bottom. The basket avoids damages at the instruments and at the tank bottom. Do not staple instruments, overloading reduces the cleaning efficiency. Pliers and scissors must be opened completely or detached, if necessary. Instruments must be completely covered with cleaning liquid, air has to escape from cavities and hoses.
- When disinfecting and cleaning contaminated medical instruments in ultrasonic baths, the hygienic safety is important after the treatment. In case of improper and irregular disinfection and cleaning a micro-biological contamination is possible caused by settlement of microorganism, especially at the tank rim or in the drain section. Cross infection might be a result. Therefore, the tank and housing surfaces must be regularly disinfected and cleaned.

According to the hygienic plan the user must carry out regular disinfection and cleaning with a surface disinfecting agent certified by DGHM*, or another efficiently classified agent.

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6. Error findings ?

- Oscillations are too weak and non-uniform, the oscillation noise is too loud; the cleaning effect is insufficient:
 - Has the liquid been properly degassed? ⇒ sonicate for 15 minutes.
 - Is the unit overloaded? ⇒ remove some parts.
 - Non-uniform noise (wobbling) is not a defect ⇒ change the fluid level.
 - Heating defect:
No objection to operate the unit without the heating.
 - Minor erosion marks at the tank bottom? ⇒ common wear and tear effect.
The unit is alright.
-
- The unit is equipped with an excess voltage arrestor.
It protects the unit from voltage surges ⇒ mains fuse can be tripped.
 - Repairs should only be carried out by authorised service.
 - Replace defective parts only by original parts.

7. Technical data

compact units are RFI-proofed and marked !.
Compliance with the limits according to EN 55011 and EN 50082-1.

Nominal voltage: 230V~ 50/60 Hz (115 V on request)
High frequency: 35 kHz; TK 100 /H 43 kHz

Environmental conditions in accordance with EN 61 010-1 (IEC 1010-1)

Degree of soiling:	2 according to IEC 664
Excess voltage category:	II
Protection class:	IP 32
permissible ambient temperature:	5 - 40 °C
permissible relative humidity up to 31 °C:	80 %
permissible relative humidity up to 40 °C:	50 %

Dew not permitted.

Information for the use as a medical product

Denotation:	Ultrasonic cleaning device
UMDNS-nomenclature (ECRI):	14-263
Purpose determination:	Ultrasonic-disinfection and cleaning of medical and dental instruments
Classification (Medical Devices Directive 93/42/EWG, appendix IX):	Class I; active, not invasive, not implantable medical device
Type, model, serial number, year of manufacture:	see type designation plate on rear of device

Information according to IEC 601-1:

Protection class: I
Protection type (B, BF, CF): not applicable, no application part available

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UR 3 | 70.791.0005/6

EU DECLARATION OF CONFORMITY

We, represented by the undersigned, hereby declare that the above device complies with the following directives and harmonised standards:

Low Voltage Directive 2014/35/EU

Electromagnetic compatibility 2014/30/EU

Restriction of hazardous substances (RoHS) 2011/65/EU

Authorised person for compilation of the technical documentation:

Julia Kürten (Technical Documentation)

Furthermore, we declare that the relevant technical documentation for the above device has been prepared in accordance with Annex VII Part A of the Machinery Directive and we undertake to submit the documentation to the market surveillance authorities on request.

In the event of a modification of the device not agreed on by Retsch GmbH, as well as the use of non-approved spare parts or accessories, this declaration loses its validity.

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Copyright

© Copyright by
Retsch GmbH
Haan, Retsch-Allee 1-5
D-42781 Haan
Federal Republic of Germany