# CoaguChek XS PT Test

<span>[</span> REF <span>]</span>	<span><span>∇</span></span>	<span>[</span> SYSTEM <span>]</span>
<b>04625374</b>	6	CoaguChek® XS
<b>04625358</b>	24	CoaguChek® XS Plus
<b>04625315</b>	2x24	CoaguChek® XS Pro

**English**

**Intended use**

Test strips for the quantitative in vitro determination of prothrombin time<sup>(a)</sup> in capillary blood or from non-anticoagulated venous blood using CoaguChek XS/XS Plus/XS Pro meters.

Suitable for self-testing.

a) Also known as Quick’s test or Quick’s value or under the abbreviation PT

**Test principle**

Electrochemical measurement of prothrombin time following activation of blood coagulation with human recombinant thromboplastin. Each test strip has a test area containing a prothrombin reagent. When blood is applied, the reagent is solved and an electrochemical reaction takes place which is transformed into a clotting time value being displayed on the meter screen.<sup>1,2</sup>

These instructions feature two symbols that draw your attention to important information.

⚠ This symbol indicates that the result could possibly be incorrect or that you could be at risk of damaging your health.

ℹ This symbol draws your attention to other important information.

**Measuring range and therapeutic range**

The CoaguChek XS/XS Plus/XS Pro meter can display results in the following units: INR (International Normalized Ratio), % Quick and Seconds.

The measuring ranges refer to the system and are defined by the technical limits of the meter and test strip. The measuring ranges are as follows: INR: 0.8-8.0 / % Quick: 120-5 / Seconds: 9.6-96

Your treating physician will tell you what your individual therapeutic range is.

The relationship of INR to % Quick may vary slightly depending on the particular lot of CoaguChek XS PT Test.

⚠ **How to proceed if test results are outside the therapeutic range**

If the measured PT result is unusually high or low repeat the test. If the PT result is still outside the therapeutic range specified by your treating physician, immediately contact your physician and ask for the appropriate (anticoagulant) measures to take in order to reduce risks that could be encountered due to excessive anticoagulation (danger of bleeding) or insufficient anticoagulation (risk of thrombosis).

**Notes on PT self-testing**

In cases where PT self-testing by the patient supplements physician care, the following should be noted: PT self-testing supplements physician care, but is not a substitute for it. PT self-testing gives the anticoagulated patient more security in everyday life. Results should be recorded in a notebook (patient’s logbook) and presented to the physician at every consultation. This makes it easier for the physician to assess the overall quality of therapy management.

**Storage and stability**

Store at 2-30 °C.

The test strips can be used up until the expiry date printed on the box and test strip container.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

⚠ This is required in order to prevent the remaining test strips from deteriorating through exposure to external influences like humidity.

**Materials provided**

- Test strips and 1 code chip

**Materials required (but not provided)**

- [REF] 04800842, CoaguChek XS Plus meter or
- [REF] 05530199, CoaguChek XS Pro meter or
- [REF] 04625412, CoaguChek XS meter
- [REF] 04696522, CoaguChek XS PT Controls
- Lancing device, e.g. CoaguChek XS Softclix®
- Lancets (e.g. CoaguChek Softclix Lancet or Accu-Chek® Safe-T-Pro Plus)

⚠ The CoaguChek XS Softclix lancing device has been specifically developed for self-monitoring by the patient. **It is not suitable for use in hospitals or for testing different persons as there is a risk of contamination.** Health care professionals performing tests on more than one patient must be aware that there is a potential risk of infection. Any object coming into contact with human blood is a potential source of infection.<sup>3</sup>

**Sample material**

Use only fresh capillary blood or venous whole blood with no anticoagulants (heparin, EDTA, citrate, oxalate, or other substances) added. If you use capillaries, do not use glass capillary tubes or capillary tubes that contain anticoagulants.

**Before you test**

Before carrying out your first blood test, please read the meter manual carefully to acquaint yourself with the way it works.

Use the CoaguChek XS PT Test test strips with the CoaguChek XS/XS Plus/XS Pro meter only.

ℹ **Getting ready**

Each code chip belongs to a particular lot of test strips. Make sure that the number on the code chip matches the number on the label of the test strip container. If you have opened a new pack of CoaguChek XS PT Test test strips, you must replace the old code chip in the meter with the one included in the new pack.

**For use with CoaguChek XS:**

Every time you insert a test strip, the meter display shows the number of the inserted code chip. Compare this number with the code printed on the test strip container. When they are identical, confirm by pressing the M-button. Keep the code chip in the meter as long as the respective test strip lot is used.

⚠ Using the wrong code chip can produce incorrect results.

**For use with CoaguChek XS Plus/XS Pro:**

When a test strip lot is used for the first time, the meter requests the corresponding code chip. If test strip lot and code chip do not correspond, an error message is displayed and a measurement is not possible.

- Perform the test as described in the meter manual and in this method sheet. All operating steps are described and illustrated in detail in the meter manual.
- Carefully wash your hands with soap and warm water. Then thoroughly dry them with a clean towel before you lance a fingertip with the lancing device to draw the drop of blood required.
- Prepare the lancing device by inserting a **new** lancet. Prime the lancing device. **Use lancets only once** because of the possible infection risk. A used lancet is a perfect breeding site for germs. Please refer to the instructions for use that come with the lancing device and the CoaguChek XS/XS Plus/XS Pro meter.

**Collecting sample and testing**

ℹ **Before each test make sure that the correct code chip for that lot of test strips is inserted in the meter.** Place the meter on a level surface or hold it in your hand so that it is roughly horizontal. Use the test strip within 10 minutes from the time you took it out of its container.

Tightly re-cap the container immediately after removing a test strip.

**Collecting a blood sample**

Test fresh capillary or venous blood immediatly. When using capillary blood apply the sample to the test strip within 15 seconds after lancing the fingertip and within 30 seconds when using venous blood. Applying blood later than that may produce an inaccurate result as the coagulation process will have begun.

⚠ **Please note the following when using capillary blood from a fingertip**

Apply the first drop of blood to the CoaguChek XS PT Test test strip. Never add more blood to the test strip after the test has begun or perform another test with blood from the same puncture site. You can either apply the drop of blood to the sample application area of the strip from above, or hold it against the side of the sample application area. The test strip draws up the blood by capillary action. You hear a beep tone when you have applied enough blood (provided the beep tone lasted long enough). Perform the test exactly as described in the meter manual.

ℹ **Do not touch or remove the CoaguChek XS PT Test test strip when a test is in progress. Wait until the result is displayed.** Record the test result in your logbook. After testing, you can discard the used lancet and CoaguChek XS PT Test test strip with your regular household waste. Dispose of used lancets carefully (e.g. use a sturdy sharps container with lid) so that the needles cannot cause injury to yourself or to others. **Healthcare professionals:** Dispose of used items in line with the disposal and infection prevention policy of your hospital, institute or medical practice.

**Quality control**

The CoaguChek XS/XS Plus/XS Pro meters have a number of inbuilt system checks. For details, please refer to the meter manual. The test strip has an integrated quality control function. Quality control and system checks with test fluids as used in certain other systems are generally not required by the CoaguChek XS/XS Plus/XS Pro meters.

Follow the applicable government regulations and local guidelines for quality control.

⚠ **Products only deliver reliable results when they are used as intended.**

**Cleaning and disinfection instructions**

To avoid malfunction of the meter, please follow the instructions provided in the according meter manual section.

**Possible causes of errors**

If problems occur during testing, please check the following:

- Are you testing exactly as described in the meter manual for your CoaguChek XS/XS Plus/XS Pro meter and in this method sheet? Please read the instructions carefully.
- Have you correctly stored the CoaguChek XS PT Test test strips (see the "Stability and storage" section of this method sheet)? Always make sure the test strips are correctly stored and that you carry out the test within 10 minutes of removing a test strip from its container.
- Is the test strip guide dirty? Clean the meter as described in the meter manual. Repeat the test with a new test strip.
- Errors displayed on the meter are generally due to an activation of the system fail safe mechanisms which are designed to prevent the release of wrong measurement results.

Special attention is required if the following errors occur:	
<span> </span>	Only valid for serial numbers: < UPO0400000 / < UQ0040000 / < U76001453
<b>“error 4”</b> <b>“error 8”</b>	Please contact your local customer support and service center if <ul style="list-style-type: none"><li>“error 4” appears when turning on the meter first time after cleaning</li> <li>“error 8” appears during the first measurement after cleaning</li></ul>

Special attention is required if the following errors occur:	
<b>“error 6”</b> <b>“error 7”</b>	In rare cases, these errors may also be received for patients under clinical conditions which can lead to abnormal or unusually long clotting times. This can happen, for example, during treatment with vitamin K antagonists in combination with antibiotics, and/or chemotherapeutics, or extremely high concentrations of oxidizing substances, e.g. after a vitamin C infusion.
<b>"E-406"</b> XS Plus/ XS Pro from SW version 03.01.00 onward	If one of these errors appears again after the test has been repeated, contact your physician without delay to get the value checked using a laboratory method.

- If you receive other error messages, please refer either to the meter manual or contact your local customer support and service center (for contact details see your meter manual).

**Additional information for healthcare professionals**

**Intended use**

For determination of prothrombin time (PT/INR) in patients on oral anticoagulant therapy with vitamin K antagonists.<sup>4,5,6,7</sup>

**Reagents**

The test strip contains human recombinant thromboplastin, stabilizers and preservatives.

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

**Calibration**

Each lot of test strips is calibrated to a reference lot that is traceable to the WHO International Reference Preparations. For the purpose of providing universal INR results, the Mean Normal Prothrombin Time (MNPT) has been established as 12 seconds for healthy volunteers and the International Sensitivity Index (ISI) for the system has been established as 1.

**Limitations – interference**

The blood drop must be a minimum of 8 µL in volume. Low sample volume will cause an error message.

Testing performed with the following in vitro spiked samples or native blood samples (or testing of Triglycerides) indicated no significant effect on test results:

- Ascorbic acid up to 170 µmol/L (30 mg/L)
- Bilirubin up to 513 µmol/L (30 mg/dL)
- Hemolysis up to 0.62 mmol/L (1000 mg/dL)
- Triglycerides up to 5.7 mmol/L (500 mg/dL)

Hematocrit ranges between 25 % and 55 % do not significantly affect test results.

The CoaguChek XS PT Test is insensitive to unfractionated and fractionated heparin concentrations up to 2 IU/mL blood. Information about pharmacokinetics and the dosage-plasma-relation of the respective heparin are provided in the package insert by the legal manufacturer of the drug.

PT results may be influenced from patients treated with Direct Oral Anticoagulants (DOACs) e.g. rivaroxaban, apixaban, edoxaban, betrixaban and dabigatran or anticoagulants other than vitamin K antagonists (e.g. hirudin and other thrombin inhibitors). For such patients, medical decisions should not be based on CoaguChek system measurements.

**Note:** Samples from patients treated with the following drugs must not be tested with the system: protamine sulfate, oritavancin, calcium dobesilate, fondaparinux.

⚠ The action of oral anticoagulants (coumarin derivatives) can be increased or weakened when other medication is taken simultaneously (e.g. antibiotics, but also prescription-free medication like pain relievers, antirheumatic medication and medication against influenza). This, in turn, can also lead to either an increase or a decrease in prothrombin time (INR). Additional medication should only be taken if prescribed by your physician. If other medication is taken, it is recommended that the prothrombin time be checked more frequently and that the anticoagulant dose be subsequently adjusted as directed by the treating physician. Anti-phospholipid antibodies (APA) like Lupus antibodies (LA) may falsely prolong coagulation times, i.e. they may cause false-high INR values and false-low Quick values. Where APA are known to be present, it is imperative that a result be obtained using an APA-insensitive laboratory method for comparison. Hirudin is not neutralized and leads to false-high INR values and false-low Quick values.<sup>8</sup>

Drug interferences are measured based on recommendations given in CLSI guidelines EP07 and EP37 and other published literature. Effects of concentrations exceeding these recommendations have not been characterized.

For diagnostic purposes, the results should always be assessed in conjunction with the patient’s medical history, clinical examination and other findings.

**Specific performance data**

**Precision**

The precision of results using capillary or venous blood was INR CV < 4.5 % using capillary blood, and < 3.5 % using venous blood in the normal as well as in the therapeutic range.

**Method comparison**

Clinical studies were conducted in which venous and capillary blood results from the CoaguChek XS/XS Plus/XS Pro Systems were compared with venous blood results obtained using the laboratory reference method Innovin (Dade-Behring). The majority of slopes was found between 0.93 and 1.04 for venous results, and between 0.92 and 1.03 for capillary results.<sup>9,10,11</sup>

**Note:** While the INR unit was created in general for the best possible comparison between methods (Point-of-Care vs. laboratory methods and between laboratory methods), several factors can have a significant or even systematic influence on the comparison of PT/INR results obtained with different methods. The most important factor, when changing methods, is the type of thromboplastin used (i.e. human recombinant, rabbit or bovine).<sup>9</sup>

The CoaguChek method uses human recombinant thromboplastin. Therefore, the comparability to tests using other human recombinant thromboplastins is best, whereas higher deviations can occur with other thromboplastin types.

However, those higher differences between thromboplastins of different (rabbit, bovine) origin are not an issue specific for CoaguChek assays. Similar differences can be observed when a human recombinant thromboplastin-based laboratory method is compared with several other (rabbit, bovine) laboratory methods. To minimize these differences, in a monitoring situation, it is recommended that each laboratory uses results from a method using only one type of thromboplastin for each patient.

If testing methods for the comparison of patient INR values are changed, especially when methods with thromboplastins of a different origin are used, the deviations which may occur need to be taken into account.

Other potential influencing factors:

- Variations between different laboratory instruments and between different reagent lots
- Pre-analytics, i.e. sample tubes of different manufacturers used for laboratory testing

Sample collection tubes other than Sarstedt for the PT INR determination with the reference method may lead to different results in method comparison.<sup>12,13</sup>

For further information, please refer to the appropriate meter manual and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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- Leichsenring I, Plesch W, Unkrig V, et al. Multicentre ISI assignment and calibration of the INR measuring range of a new point-of-care system designed for home monitoring of oral anticoagulation therapy. Thromb Haemost 97:846-861, 2007.
- Plesch W, Van den Besselaar AMHP: Validation of the international normalized ratio (INR) in a new point-of-care system designed for home monitoring of oral anticoagulation therapy. Int Jnl Lab Hem 31:20-25, 2009.
- Clinical and Laboratory Standards Institute: Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline – Third Edition; CLSI document M29-A3, 2005.
- Ansell J, Jacobson A, Levy J et al. Guidelines for implementation of patient self-testing and patient self-management of oral anticoagulation. International consensus guidelines prepared by International Self-Monitoring Association for Oral Anticoagulation. Int J Cardiol 2005;99:37-45.


- Menéndez-Jándula B, Souto J et al. Comparing Self-Management of Oral Anticoagulant Therapy with Clinic Management. Ann Intern Med 2005;142:1-10.
- Taborski U, Wittstamm F-J, Bernardo A: Cost-effectiveness of self-managed anticoagulant therapy in Germany. Sem Thromb Hemostas 1999;25:103-109.
- Tripodi A, Breukink-Engbers WGM, van den Besselaar AMHP: Oral Anticoagulant Monitoring by Laboratory or Near-Patient Testing: What a Clinician Should Be Aware Of. Sem Vasc Med 2003;3:243-254.
- Moll S and Ortel TL Monitoring Warfarin Therapy in Patients with Lupus Anticoagulants. Annals of Internal Medicine 1997;127:177-185.
- Plesch W, Wolf T, Breitenbeck N et al. Results of the performance verification of the CoaguChek XS System. Thromb Res 2008;123:381-389.
- Deom A, Reber G, Tsakiris DA et al. Evaluation of the CoaguChek XS Plus system in a Swiss community setting. Thromb Haemost 2009;101:988-990.
- Multicenter Performance Evaluation of the CoaguChek XS System; Evaluation Report from Roche Diagnostics available upon request.
- van den Besselaar AMPH, Hoekstra MMCL, Witteveen E, et al. Influence of Blood Collection Systems on the Prothrombin Time and International Sensitivity Index Determined With Human and Rabbit Thromboplastin Reagents. Am J Ciin Pathol 127:724-729, 2007.
- van den Besselaar AMPH, Rutten WPF, Witteveen E. Effect of magnesium contamination in evacuated blood collection tubes on the prothrombin time test and ISI calibration using recombinant human thromboplastin and different types of coagulometer.Thrombosis Research 115:239-244, 2005.

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
	COAGUCHEK, ACCU-CHEK and SOFTCLIX are trademarks of Roche.
	All other product names and trademarks are the property of their respective owners.
	Additions, deletions or changes are indicated by a change bar in the margin.
	© 2020, Roche Diagnostics
	

www.coagucheck.com

	Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim www.roche.com