



06879942001V6.0

Roche CARDIAC D-Dimer

cobas[®]

REF 04877802 190

▽ 10

SYSTEM cobas h 232

English

Intended use

Quantitative immunological test for the detection of d-dimer in heparinised venous blood for use with the **cobas h 232** instrument.

The Roche CARDIAC D-Dimer test serves as an aid when deep venous thrombosis and pulmonary embolism is suspected. A negative d-dimer result is an indication that these diseases can be ruled out with high probability.

Summary

D-dimer is a degradation product of crosslinked fibrin. The d-dimer concentration is a measure of the fibrinolytic activity of plasmin in the vascular system. Elevated concentrations of d-dimer indicate increased coagulatory and fibrinolytic activity. With a normal d-dimer value, acute deep vein thrombosis and pulmonary embolisms may be ruled out with very high reliability.^{1,2,3,4,5,6,7,8,9,10,11,12,13}

Test principle

The test contains two monoclonal antibodies against fibrin degradation products which contain the d-dimer structure element. One of the antibodies is gold-labelled, the other biotinylated. The antibodies form a sandwich complex with the d-dimer in the blood. Following removal of erythrocytes from the sample, plasma passes through the detection zone in which the gold-labelled d-dimer sandwich complexes accumulate and the positive signal is displayed as a reddish line (the signal line). Excess gold-labelled antibodies accumulate along the control line, signalling that the test was valid. The intensity of the signal line increases in proportion to the d-dimer concentration.

The optical system of the instrument detects the two lines and measures the intensity of the signal line. The integrated software converts the signal intensity to a quantitative result and shows it in the display.

Reagents

One test contains:

Biotinylated mouse monoclonal anti-d-dimer antibodies $\geq 1.0 \mu\text{g}$

Gold-labelled mouse monoclonal anti-d-dimer antibodies $\geq 1.0 \mu\text{g}$

Buffer and non-reactive components $\geq 2.8 \text{ mg}$

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

Storage and stability

Until the printed expiration date at 2-8 °C.

Up to 1 week at room temperature (15-25 °C).

The test can be used immediately after removal from the refrigerator.

The test must be used within 15 minutes once the pouch has been opened.

Sample stability: 8 hours at room temperature. Do not refrigerate or freeze sample.

Specimen collection and preparation

Use **heparinised venous whole blood** only.

Do not use other anticoagulants, capillary blood, serum or plasma, blood collection tubes containing EDTA, citrate, sodium fluoride or other additives.

The following heparin blood collection tubes have been tested: Sarstedt Monovette, Becton Dickinson Vacutainer, Becton Dickinson Vacutainer PST II, Greiner Vacuette, Terumo Venosafe. In the case of Sarstedt Monovettes, only tubes without separating gel are suitable.

No data is available for blood collection tubes supplied by other manufacturers. An influence on the test result in individual cases cannot be ruled out.

Sample volume: 150 μL

Materials provided

- REF 04877802190, Roche CARDIAC D-Dimer test
- 1 code chip

Materials required (but not provided)

- REF 11622889190, Roche CARDIAC Pipettes, 20 disposable syringes (150 μL)
- REF 04890523190, Roche CARDIAC Control D-Dimer (2 x 1 mL)
- REF 04880668190, Roche CARDIAC IQC
- REF 04901126190, **cobas h 232** instrument
- REF 04901142190, **cobas h 232** instrument with scanner
- General laboratory equipment

Calibration

The Roche CARDIAC D-Dimer test is calibrated against the Tina-quant D-Dimer test using citrate plasma.

The instrument automatically reads in the lot-specific calibration data from the code chip, eliminating the need for calibration by the user.

Lot code

Every kit contains a lot-specific code chip. The instrument display prompts the user to insert the chip. To ensure that the code chip and test strip lot match, compare the lot number in the display with the number on the code chip. The code chip provides the instrument with all required lot-specific information. An error message is displayed if the wrong code chip is inserted for a test strip lot.

Quality control

For quality control, use Roche CARDIAC Control D-Dimer.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The instrument automatically calculates the concentration of each sample.

The reaction time for the Roche CARDIAC D-Dimer test to display a quantitative result is 8 minutes. In addition, approximately 2 minutes are required for sample detection.

Limitations - interference

The assay is unaffected by icterus (bilirubin $\leq 20 \text{ mg/dL}$), hemolysis (Hb $\leq 200 \text{ mg/dL}$), lipemia (triglycerides $\leq 470 \text{ mg/dL}$), haematocrit values in the range of 26-56 %, and biotin $\leq 30 \text{ ng/mL}$.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. $> 5 \text{ mg/day}$) until at least 8 hours following the last biotin administration.

No interference was observed from rheumatoid factors up to a concentration of 300 IU/mL.

High concentrations of lipoic acid (e. g. in pharmaceuticals or as food additive) can lead to lower measurement values.

Very high d-dimer concentrations (approx. $> 50 \mu\text{g/mL}$) can lead to lowered values (hook effect). However, this effect does not result in false normal values. Alternatively, the control line may fail to appear and the instrument may display an error message. In this case, determination must be carried out using another method, like the Tina-quant D-Dimer test from Roche.

High concentrations of d-fragments as may appear under thrombolytic therapy can lead to lower measurement values.

Patient samples may contain heterophilic antibodies which could react in immunoassays to give falsely elevated or decreased results. Reasons for the presence of heterophilic antibodies might be for example elevated levels of rheumatoid factors or the treatment of patients with monoclonal mouse antibodies for therapeutic or diagnostic purposes.

The Roche CARDIAC D-Dimer test contains ingredients that minimise interference from heterophilic antibodies. However, complete elimination of interference from all samples cannot be guaranteed. Interferences caused by pharmaceuticals at therapeutic concentrations are not known.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.



Measuring range

0.1-4 µg/mL.

1 µg/mL corresponds to 1 µg FEU/mL.

Expected values

The normal range for the Roche CARDIAC D-Dimer test includes values up to 0.5 µg/mL. Values above 0.5 µg/mL are to be considered as pathologically elevated.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the instruments are given below. Results obtained in individual laboratories may differ.

Precision

Repeatability was measured with 3 lots of Roche CARDIAC D-Dimer tests and heparinised human blood. The majority of the variation coefficients was below 11 %. Intermediate precision was measured with the Roche CARDIAC Control D-Dimer quality control in 5 different hospitals. The majority of the variation coefficients was below 10 % (level 1) and below 17 % (level 2).

Method comparison

A comparison using 3 different lots of the Roche CARDIAC D-Dimer test with the Tina-quant D-Dimer test in a clinical patient population showed slopes between 0.94 and 1.03 in the majority of the method comparisons. The majority of the correlations in these method comparisons were ≥ 0.93 .

References

- 1 Baker WF. Diagnosis of deep venous thrombosis and pulmonary embolism. *Med Clin North Am* 1998;82:459-476.
- 2 Becker DM, Philbrick JT, Bachhuber TL, et al. D-dimer testing and acute venous thromboembolism. A shortcut to accurate diagnosis? *Arch Intern Med* 1996;156:939-946.
- 3 Bernardi E, Prandoni P, Lensing AW, et al. on the behalf of the Multicentre Italian D-dimer Ultrasound Study Investigator Group. D-dimer testing as an adjunct to ultrasonography in patients with clinically suspected deep vein thrombosis: prospective cohort study. *BMJ* 1998;317:1037-1040.
- 4 Bounameaux H, de Moerloose P, Perrier A, et al. Plasma measurement of D-dimer as diagnostic aid in suspected venous thromboembolism: an overview. *Thromb Haemost* 1994;71:1-6.
- 5 Brill-Edwards P, Lee A. D-dimer testing in the diagnosis of acute venous thromboembolism. *Thromb Haemost* 1999;82:688-694.
- 6 Demarmels Biasutti F, Lämmle B. Beitrag des Hämostaselabors bei der Diagnostik der tiefen Venenthrombose. *Therap Umschau* 1996;53:265-271.
- 7 Dempfle CE, Hafner G, Lestin HG, et al. Multizentrische Evaluierung von Tina-quant [a] D-Dimer. *J Lab Med* 1996;20:31-37.
- 8 Dempfle CE, Schraml M, Besenthal I, et al. Multicentre evaluation of a new point-of-care test for the quantitative determination of D-dimer. *Chim Clin Acta*.
- 9 van der Graaf F, van den Borne H, van der Kolk M, et al. Exclusion of deep vein thrombosis with d-dimer testing – Comparison of 13 d-dimer methods in 99 outpatients suspected of deep venous thrombosis using venography as reference standard. *Thromb Haemost* 2000;83:191-198.
- 10 Knecht MF. Stellenwert der D-Dimere in der Diagnostik thromboembolischer Erkrankungen. *medwelt* 1997;48:52-65.
- 11 Lestin HG, Ehrenteich U, Hergert M, et al. Untersuchungen zur Relevanz hämostaseeologischer Marker für die Diagnostik, Therapie und Verlaufskontrolle tiefer Venenthrombosen. *Clin Lab* 1996;42:745-756.
- 12 de Moerloose P, Michiels JJ, Bounameaux H. The place of D-Dimer testing in an integrated approach of patients suspected of pulmonary embolism. *Semin Thromb Hemostasis* 1998;24:409-412.
- 13 Perrier A, Buswell L, Bounameaux H, et al. Cost-effectiveness of noninvasive diagnostic aids in suspected pulmonary embolism. *Arch Intern Med* 1997;157:2309-2316.

For further information, please refer to the appropriate operator's manual for the instrument concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

SYSTEM

Analyzers/Instruments on which reagents can be used

COBAS, COBAS H, IQC, ROCHE CARDIAC and TINA-QUANT are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2014, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

