WESTMED

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Please read this manual carefully before using the device. Failure to follow these instructions could result in incorrect measurements or damage the device.

Version of the manual: Ver 1.1

Last update: 2021/06/30

All rights reserved.

### **Please note**

- The contents of this manual are subject to change without notice.
- The information provided by the company is believed to be accurate and reliable. However, the company assumes no responsibility for the use of the manual, or any infringements of patients' or other rights of third parties that may result from its use.

### **Instructions for safe operation**

- Check the device to ensure there is no visible damage that may affect the user's safety or the measurement performance of the sensors

and clips. The company recommends inspecting the device before each use. If there is any obvious damage, do not use the device.

- Particular caution is advised if the Oximeter is used constantly at ambient temperatures over 37°C, as it could cause burns if the sensor becomes overheated.
- Any maintenance necessary may only be performed by qualified service technicians. Users are not permitted to service this device.
- The Oximeter may not be used together with devices and accessories not specified in this manual.

### Safety precautions

- Danger of explosion — DO NOT use the Oximeter in environments that could contain inflammable gases, such as certain types of inflammable anesthetics.
- DO NOT use the Oximeter while the patient is underground MRI or CT scanning. This device is NOT MRI-compatible.

## Warnings

- Discomfort or pain may occur if the Oximeter is used continuously on the same part of the body for a prolonged period, especially if the patient suffers from poor microcirculation. The company recommends not applying that the Oximeter to the same physical location ((oder)) part of the body for longer than 2 hours. If any abnormality is noticed, please change the position of the Oximeter.
- DO NOT attach this device to swollen or tender tissue.
- Infrared light is invisible. The light emitted from the device is harmful to the eyes. Do not stare directly into the light.
- The Oximeter is not intended for use as a treatment device.
- Local laws and regulations must be followed when disposing of the device.

## Safety information

- Keep the Oximeter away from dust, vibration, corrosive substances, explosive materials, high temperatures, and moisture.

- The device should be kept out of the reach of children.
- If the Oximeter becomes wet, please stop using it and do not resume operation until it is dry and has been checked for correct operation. If the device is taken from a cold environment to a warm, humid environment, please do not use it immediately. Allow at least 15 minutes for the Oximeter to reach the ambient temperature.
- DO NOT press any of the buttons on the front panel with sharp or pointed objects.
- DO NOT disinfect the Oximeter using high temperatures or high-pressure steam.
- See section 7 for instructions regarding cleaning and disinfection.
- The device complies with protection class IP22 and is protected against the intrusion of harmful solid matter and liquids. The device is therefore protected against solid foreign matter 12.5 mm or greater in size and against vertically falling water drops when the casing is tilted by up to 15°.

- Please protect the device from the effects of lint, dust, light (including sunlight), etc.

### **Declaration of Conformity**

The manufacturer hereby declares that this device complies with the following standards:

IEC 60601-1: 2005 + A1: 2012 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.

ISO 80601-2-61:2017 – Medical electrical equipment – Part 2-61: Particular requirements for the basic safety and essential performance of pulse oximeter equipment.

And it also follows the provisions of the council directive MDD 93/42/EEC.

**Caution: U.S. federal law restricts the sale or use of this device by or on the order of a physician.**

## Intended use

This fingertip oximeter is intended for measuring the pulse rate and functional oxygen saturation (SpO<sub>2</sub>) via a patient's finger. It is applicable for randomly checking the oxygen saturation and pulse rate of adult and child patients at home and in medical clinics. Models with external sensors and threshold indicators may be used for longer periods of time, depending on the suitability of the sensor selected.

## Description

- **Muting the sound**  
For models fitted with both threshold indicators and pulse beep functions, when both the beep and the threshold indicator signal are activated, the Display key will then also function as the Mute key. Press the Display key for a short time to switch off the threshold indicator signal and the pulse beep for 90 seconds.
- **Measuring mode** 1.  
**Random checking mode:** The measurement

begins automatically when the finger is correctly inserted into the finger clip. The measuring time lasts 30 seconds and is accompanied by a countdown timer. After 30 seconds, the readings for the SpO<sub>2</sub>, pulse, and pulse rhythm analysis will be displayed on the screen. Once the finger is removed from the clip, the display will be cleared and the Oximeter shuts down automatically.

2. **Continuous mode:** The measurement begins automatically when the finger is correctly inserted into the finger. The causing the Oximeter to shut down automatically.

- **Record list**

1. When the Oximeter is switched off, it automatically records a single group of stable readings in the record list, regardless of the measurement mode in which they were recorded. However, if the time from displaying the current readings to the end of measurement is less than 5 seconds, no reading will be stored.

2. Up to 12 groups of records can be stored in the record list, the newest record as M1, and

the oldest record is marked as M12. The new record will override the previous record.

3. If the batteries are removed from the device, the records will be deleted.
4. When the device is switched off, hold down the Display key to view the record recall screen. When the record recall screen is displayed, briefly press the Display key to scroll through the records. If you do not touch a key for 6 seconds, the Oximeter will automatically switch off.

### **Battery installation**

1. Insert two AAA batteries into the battery compartment in accordance with the polarity markings.


2. Replace the cover.

Please ensure that the batteries are correctly installed, otherwise the device will not work.

Make sure to remove batteries if the device is not being used for more than 7 days to prevent any potential damage caused by battery leakage. Any such damage is not covered under the product warranty.

## Operation

Open the clip and place the finger between the rubber cushions of the clip making sure the finger is in the correct position), and then close the clip over the finger.

- After 2 seconds, the Oximeter will switch on automatically and start to measure
- The data display screen will then be shown

## Note

When “- -” displays on device screen, it indicates the reading is unavailable right now. It may be caused by:

- Excessive movement;
- Poor signal;
- Ambient light;

Usually, the reading will recover in a few seconds when at rest.

## Display

Briefly press the Display key to change the direction of the display.

## Menu

Press and hold the Display key to open the settings menu. You can scroll through the menu by briefly pressing the Display key. To change a setting, press and hold the display key at the relevant menu item. The setting to be changed should now flash. The setting can now be changed by briefly pressing the Display key. To confirm the change, press and hold the Display key again. To activate the changes and exit the menu, scroll to the menu option “Save exit menu” and keep the Display key pressed.

## Information for correct use

- The finger should be pcorrectly inserted into the sensor.
- The finger may not move and must be relaxed during measurement.
- Do not insert a wet finger directly into the sensor.
- Avoid placing the device on the same limb

that is wrapped with a cuff for blood pressure measurement or during venous infusion.

- Ensure that the light emitted by the device is not blocked in any way, e.g. the fingernail may not be varnished or painted.
- Vigorous exercise and electrosurgical device interference may affect the measuring accuracy.
- The orientation sensor works on the basis of the earth's gravity. A small movable metal ball is built in the orientation sensor in order to detect the orientation of the Oximeter. If you wish to change the Oximeter's display direction, if you move the Oximeter too slowly, the movable metal ball will also move slowly due to too little acceleration. Consequently, the orientation detection response will be delayed. Acceleration is required in order for the orientation sensor to sense the orientation change with sufficient speed.
- Nail polish may affect the measuring accuracy. If the fingernail is too long, it may cause the measurement to fail or deliver an inaccurate result.

- Any highly intensive lighting, red lights, infrared heaters, or strong sunshine, etc. may cause inaccurate measurement results. Please put an opaque cover on the sensor or perform the measurement at a different location if necessary.
- If the first reading produces a poor waveform (irregular or not smooth), then the reading is most likely incorrect. You will receive a more stable reading if you wait for a moment or restart the device if necessary.

## Technical specifications

### 1. SpO<sub>2</sub> measurement

Transducer: dual-wavelength LED sensor with the following wavelength:

Red light: 663 nm, infrared light: 890 nm.

Maximum average optical output power:

≤2 mW

SpO<sub>2</sub> display range: 35%–100%

SpO<sub>2</sub> measuring accuracy:

≤2% for SpO<sub>2</sub> range from 70% to 100%

## 2. Pulse rate measurement

PR display range: 30 bpm~240 bpm

PR measuring accuracy:  $\pm 2$  bpm or  $\pm 2\%$

(whichever is greater)

## 3. Perfusion Index (PI) display range

0%–20%

## 4. Preset upper limits

SpO<sub>2</sub> lower limit: 90%

Pulse rate:                   Upper limit: 120 bpm  
  Lower limit: 50 bpm

## 5. Power supply requirement

2 x LR03 (AAA) alkaline batteries

Supply voltage: 3.0

Operating current:  $\leq 40$  mA

## 6. Environmental conditions:

Operating temperature:           5°C–40°C

Operating humidity:               30%–80%

Atmospheric pressure:           70 kPa–106  
kPa

## 7. Low perfusion performance:

The accuracy of SpO<sub>2</sub> and PR measurements still meets the above-mentioned precision when the modulation amplitude is as low as 0.6%.

## 8. Ambient light interference:

The difference between the SpO<sub>2</sub> values measured in natural indoor light conditions and those in darkness is less than  $\pm 1\%$ .

## 9. Dimensions:

59 mm (L) × 34 mm (W) × 30 mm (H)

Net weight: approx. 60 g

## 10. Classification

Type of protection against electric shock:

Internally powered equipment.

Degree of protection against electric shock:

Type BF for applied parts.

Degree of protection against harmful intrusion of solid foreign matter and liquids:

This device complies with IP22 standards and is protected against the harmful intrusion of solid foreign matter and liquids.

Electromagnetic compatibility: Group I, Class B.

## Scope of delivery

- Fingertip Oximeter
- User manual
- Batteries
- Lanyard

## Repair and maintenance

The expected service life (not a warranty) of this device is 5 years. In order to ensure its long service life, the device must be properly maintained.

- Please change the batteries when the low-indicator lights up.
- Please clean the surface of the device with 75% alcohol wipes before use, then let it air-dry or wipe it dry. Do not allow liquid to enter the device.
- Please remove the batteries if the Oximeter will not be used for longer than 7 days.
- Recommended storage conditions:  
Ambient temperature:  $-20^{\circ}\text{C}$  to approx.  $60^{\circ}\text{C}$ , relative humidity 10% to approx. 95%, atmospheric pressure: 50 kPa to approx. 107.4 kPa.

- The Oximeter is calibrated in the factory before sale, so there is no need to calibrate it during its life cycle. SpO<sub>2</sub> simulators may not be used to validate the accuracy of the Oximeter, but only as functional testers to verify its precision. The SpO<sub>2</sub> accuracy claimed in this manual is supported by the clinical study conducted by inducing hypoxia on healthy, non-smoking, light- to dark-skinned subjects in an independent research laboratory.
- If it is necessary to verify the precision of the Oximeter on a regular basis, the user can do so by means of an SpO<sub>2</sub> simulator, or it can be done by a locally based testing institution. Please note that the specific calibration curve (so-called R-curve) should be selected when using SpO<sub>2</sub> simulators, e.g. for the Index 2 series SpO<sub>2</sub> simulator from Fluke Biomedical Corporation, please set “Make” to “DownLoadMake: KRK,” then the user can use this particular R-curve to test the Oximeter. If the SpO<sub>2</sub> simulator does not provide the “KRK” R-curve, please ask the manufacturer for

helping to download the given R-curve to the SpO<sub>2</sub> simulator.

- **High-pressure sterilization may not be used on the device.**
- **Do not immerse the device in liquid.**
- **It is recommended that the device should be kept in a dry environment. Humidity may reduce the life of the device, or even damage it.**

### **Instructions for cleaning and disinfecting**

- Clean the surface of the sensor with a soft cloth moistened with a 75% isopropanol solution. For gentle disinfection, please use a mild bleach solution.
- Then clean the surface of the device with a cloth moistened ONLY with water and dry with a clean, soft cloth.

#### **Caution:**

Do not sterilize with a steam jet or ethylene oxide.

Do not use the Oximeter if it is damaged.

## Troubleshooting

The SpO<sub>2</sub> and pulse rate display are instable  
Position the finger correctly in the device and try again.

Device cannot be switched on

Change the batteries.

No display

Ensure that the patient is relaxed.

Direction of display does not change or is insensitive to changes

Please shake the Oximeter with a certain force to make the movable metal ball move freely. If the problem persists, the orientation sensor may not be working properly.

If the above problem persists, please contact your local service center.

## Frequently Asked Questions

Q: What is SpO<sub>2</sub>?

A: SpO<sub>2</sub> is the percentage of oxygen saturation in the blood.

Q: What is the normal range of SpO<sub>2</sub> readings for healthy people?

A: The normal range varies by individual, but is usually over 95%, otherwise, please consult your physician.

Q: What is the normal range of PR readings for healthy people?

A: Usually, the normal range is 60 bpm to approx. 100 bpm.

Q: Why do the displayed values of SpO<sub>2</sub> and PR vary over time?

A: The measured SpO<sub>2</sub> and PR values change in correspondence with changes in the patient's physiological condition.

Q: What should be done if there are no SpO<sub>2</sub> and PR readings?

A: Do not shake the finger and stay relaxed during the measurement. In addition, please avoid using the Oximeter and the cuff on the same limb used for blood pressure and oxygen saturation measurement simultaneously.

Q: How do I know that the SpO<sub>2</sub> reading is accurate?

A: Hold your breath for a while (50 seconds or longer). If the SpO<sub>2</sub> value significantly decreases, it means that the SpO<sub>2</sub> reading truly reflects the change in your physiological condition.

**Q: When do the batteries need to be replaced?**

A: The “low battery” icon will appear on the screen when the battery voltages are low. The batteries then need to be replaced.

**Q: What should I do if the Oximeter is moist or has been sprayed with water?**

A: Remove the batteries immediately and dry the Oximeter completely with a hair dryer.


**Q: What factors will affect the SpO<sub>2</sub> accuracy?**

- A: a) Intravascular dyes such as indocyanine green or methylene blue;
- b) Exposure to excessive illumination, such as surgical lamps, bilirubin lamps, fluorescent lights, infrared heating lamps, or direct sunlight;
- c) Vascular dyes or externally used products such as nail enamel or skin care products;
- d) Excessive patient movement;
- e) Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line;
- f) Exposure to high-pressure oxygen;
- g) Arterial occlusion near the sensor;

- h) Blood vessel contraction caused by peripheral vessel hyperkinesis or a drop in body temperature;
- i) Low perfusion (Perfusion Index is low).

**Please contact the local distributor or manufacturer if necessary.**

## Symbols

Symbol	Description
	CE mark
SN	Serial number
	Date of manufacture
	Attention – refer to User Manual
	UK Responsible Person
	Catalogue number
	Temperature limit
	Authorised representative in the European community
	Manufacturer (including address)
	BF type applied part
	Follow WEEE regulations for disposal
--	Unavailable readings
	Medical Device

## Appendix Analysis of results

No.	Description
1	No irregularity found
2	Suspected slightly fast pulse
3	Suspected fast pulse
4	Suspected fast pulse for short periods
5	Suspected slightly slow pulse
6	Suspected slow pulse
7	Suspected occasional short pulse intervals
8	Suspected irregular pulse intervals
9	Suspected fast pulse with short pulse intervals
10	Suspected slow pulse with short pulse intervals
11	Suspected slow pulse with irregular pulse intervals
12	Poor signal, repeat measurement

## Appendix II EMV

The device meets the requirements of IEC 60601-1-2:2014.

### **Guidance and manufacturer's declaration – electromagnetic emissions**

The Fingertip Oximeter is intended for use within the electromagnetic environment specified below. The customer or the user of the Fingertip Oximeter should assure that it is used within this environment.

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The Fingertip Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The Fingertip Oximeter is suitable for use in all facilities, including private homes and those that supply buildings with electricity for domestic purposes.
Harmonic emissions IEC61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC61000-3-3	N/A	

## Guidance and manufacturer's declaration – electromagnetic emissions

The Fingertip Oximeter is intended for use within the electromagnetic environment specified below. The customer or the user of the Fingertip Oximeter should assure that it is used within this environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment and guidance
Electrostatic discharge (ESD) IEC61000-4-2	$\pm 8$ kV contact $\pm 15$ kV air	$\pm 8$ kV contact $\pm 15$ kV air	Floors should be made of wood, concrete, or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Electrical fast transient/ burst IEC61000-4-4	$\pm 2\text{kV}$ for power Supply lines $\pm 1\text{ kV}$ input/output lines	N/A	N/A
Surge IEC61000-4-5	$\pm 1\text{kV}$ line (s) to line(s) $\pm 2\text{kV}$ line(s) to earth	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	$< 5\% U_T$ ( $> 95\%$ dip in $U_T$ ) for 0.5 cycle $< 40\% U_T$ ( $60\%$ dip in $U_T$ ) for 5 cycles $< 70\% U_T$ ( $30\%$ dip in $U_T$ ) for 25 cycles $< 5\% U_T$ ( $> 95\%$ dip in $U_T$ ) for 5 s	N/A	N/A

Power frequency (50 Hz/60 Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: $U_T$ is the AC mains voltage prior to application of the test level.			

Conducted RF IEC61000- 4-6	3 VRMS 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the Fingertip Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
-------------------------------------	--------------------------------	-----	---

<p>Conducted HF IEC61000- 4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	<p><b>Recommended separation distance</b>  <math>d = 1,2 \sqrt{P}</math>  <math>d = 1,2 \sqrt{P}</math>  80 MHz  to 800MHz  <math>d = 2,3 \sqrt{P}</math>  800MHz  to 2,5 GHz  Where P is  the maximum  output power  rating of the  transmitter  in watts (W)  according to  the transmitter  manufacturer  and d is the  recommended  separation  distance in</p>
---	--	--------------	---

			<p>meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol. </p>
--	--	--	---

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts, and TV broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Fingertip Oximeter is used exceeds the applicable RF compliance level above, the Fingertip Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Fingertip Oximeter.

b: Over the frequency range 150 kHz to 80 MHz, field strengths should be lower than 3 V/m.

Rated maximum output power of transmitter W (Watts)	Separation distance according to frequency of transmitter M (meters)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	80 MHz to 2,5 GHz
	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23

No liability is accepted for translation errors.

For transmitters rated at a maximum output power not listed above, the recommended separation distance ( $d$ ) in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## Quality Certificate

Name: Fingertip Oximeter

Model:

Date:

QA:

This product has been inspected in accordance with the standards specified in the User Manual.  
Shenzhen Creative Industry Co., Ltd.

No liability is accepted for translation errors.

Shenzhen Creative Industry Co., Ltd.  
Floor 5, BLD 9, Baiwangxin High-Tech  
Industrial Park, Songbai Road, Xili Street,  
Nanshan District, 518110 Shenzhen,  
P. R. China

Tel.: +86-755-2643 3514

Fax: +86-755-2643 0930

E-Mail: [info@creative-sz.com](mailto:info@creative-sz.com)

Webseite: [www.creative-sz.com](http://www.creative-sz.com)

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Deutschland