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ultraMEDIC GmbH & Co. KG

Im Bruch 13
56567 Neuwied · Germany

Telefon:
+49(0) 26 31/9 69 83 -0

Telefax:
+49(0) 26 31/9 69 83 -33

email:
info@ultramedic.de

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4.5. CLEANING

Cleaning can be executed applying the following indications:

METALLIC PARTS

Simply clean with a sponge drenched with water and normal neuter soap, or, if necessary, with proper disinfectant. Rinse and dry carefully. After that apply a thin protective wax coating.

PLASTIC PARTS


Clean with a sponge drenched with water diluted with disinfectant and then rinse and dry completely.

In any case, avoid products with aggressive agent (ammonia, bleach,...) and avoid utilization of abrasives to remove spots such as brushes, steel wool, knives' blades.

4.6. LUBRICATION

Lubrication can be simply executed by applying in the indicated points of the figure beneath one or two drops of medium viscosity oil.


4.4. STORING AFTER USE

The scoop stretcher, being build of rust-resisting and weatherproof materials does not require particular care if used in ordinary conditions. It is however recommended to keep it in a close place, paying attention to protect locking devices and sliding parts from dust and dirt, to obtain the best efficiency at the moment of utilization.

The scoop can be folded to reduce the required space and making storing easier.

- lay down the closed stretcher in horizontal position (A);
- place yourself beside the stretcher, near the feet zone profile
- open the right and left locking devices and unthread the two feet profiles at the same time until the folding joints are completely got out (B);
- taking hold of the frame at feet extremity turning all this part until complete superposition with the bust zone profile (C);
- store the stretcher at a proper place.

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
Stretcher opened in „A“ position

- place the top of the stretcher, already adjusted, near the head of the patient who has to be transported in such a way that the head is inside the stretcher
- insert the top of the stretcher under the head of the patient, applying the proper medical techniques of emergency
- insert the other parts of the stretcher under the patient approaching two sections until pairing and relative intervention of the locking device (feet side hook)
- secure the patient to the stretcher by means of the safety belts, adjusting their tightness and if necessary apply a cervical collar

After closing the two parts make sure about the right intervention of all locking devices (head side and feet side hook, locking of feet side profiles' extraction and safety belts)

Non use of safety belts can cause damages to the patient. Check the tightness of the belts on the patient before each utilization

The patient's transport requires the employment of almost two trained and qualified operators. In particular situations the transport with more than two operators is recommended (see point 3.3).

If you want to divide it in the two parts (pos. 1); proceed as following:

- lay the close stretcher in horizontal position;
- place yourself at the head side of the stretcher;
- holding at the same time the right and the left section, push the central button of unlocking device and separate the two sections until reaching their complete opening.
- Repeat that to unlock the feet side.

If you want to open the stretcher in „A“ position (pos. 2), unlock only the feet side hook.

4.3 PLACEMENT OF THE PATIENT AND CLOSING OF THE STRETCHER


Follow the reported procedure depending on the type of opening the stretcher.

During the placement of the stretcher pay attention to avoid to tread on the patient and/or to entangle in the stretcher clothes or hair.

Stretcher opened in two parts:

- Place the two sections, already adjusted, at the two sides of the patient who has to be transported
- Insert the two sections under the patient, applying the proper techniques of emergency, approaching the two sections until pairing and relative intervention of the locking devices (head side and feet side hook)
- Secure the patient to the stretcher by means of the safety belts, adjusting their tightness and if necessary apply a cervical collar.

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4. USE INSTRUCTIONS

The stretcher's use concerns the preliminary procedures to adjust and open it then its placement under the patient.

The minimum weight and the ergonomics allow its easy and safe shifting; the folding system permits storing with a minimum of required space.



If there is a suspicion of spinal, cervical or lumbar traumas for the patient, it is recommended an immediate removal of the patient on a proper equipment, such as spinal stretcher, removing, if necessary, the stretcher.

4.1. ADJUSTMENT OF THE STRETCHER'S LENGTH

The stretcher's length can be adjusted on 4 different positions, depending on the patient's height. Proceed as the following:

- lay the stretcher in horizontal position
- place yourself beside the stretcher, near the feet zone profile
- open left and right locking devices;
- unthread the two feet profiles and adjust in the desired position, considering that the stretcher's length can be adjusted from 1670 to 2010 mm. in 4 intervals
- once obtained the desired length close the locking devices previously opened.

Once executed the adjustment, make sure of the correct closing of locking devices checking that the feet profiles can't be unthreaded

4.2. OPENING OF THE STRETCHER

To place the stretcher under the patient without removing him, it is necessary to open the stretcher.

There are several application techniques of the scoop stretcher for example inserting one section at a time under the patient, or opening the stretcher only by feet side creating a „A“ and placing the top under the patient's head. Those applications must be chosen depending on the type of intervention.

3.3. INSTRUCTIONS FOR THE OPERATORS

The product has been made for a professional use. The operators, besides the necessary technical knowledge, must have experience specifically in the sector of assistance, particularly have reached ability in the use of the stretcher.

To work in safe and efficient conditions the stretcher requires the presence of at least two operators trained and in fit physical conditions, because they must avail themselves of such strength and muscular coordination to guarantee the maximum level of security for the patient.

In particular situations such as for example patient of excessive weight, difficult ways, presence of stair-cases, it may be necessary the employment of more than two operators. They also must verify, before and during the transport that the patient is firmly positioned in the stretcher with the safety belts.

The schedule shows the recommended number of operators for the transport relatively to patient's weight:

<u>NUMBER/KIND OF OPERATORS</u>	<u>PATIENT'S WEIGHT</u>
2 trained and qualified operators	up to 80 kgs
2 trained and qualified operators + 1 auxiliary operator	from 80 kgs to 120 kgs
2 trained and qualified operators + 2 auxiliary operator	from 120 kgs to 170 kgs
Use other devices/ vehicules for the transport	over 170 kgs

All the operators must be turned towards the patient.
In those situations which require the presence of auxiliary staff, the stretcher's control and the transport's coordination are a duty of the two qualified operators. All the employed in the transport must stay turned towards the patient. The auxiliary staff must help in lifting when necessary, following the responsible operators instructions.

In any case at least one qualified operator must never let the patient alone during all the time he stays on the stretcher.

To break the safety conditions here reported may result dangerous.

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
3. SAFETY

3.1. DIRECTIVES AND REGULATIONS

The stretcher has been projected and built according to Directive 93/42/CEE concerning medical devices, the scoop stretcher is considered class I device.

Conformity with the safety requirements provided for by Directive 93/42/EEC was obtained by applying the harmonized standard DIN EN 1865.

3.2. PROPER AND IMPROPER USE

The product is a professional one and it is intended for an employment in a medical or private situation, to make possible the transport in horizontal position of patients temporarily or permanently unable to move, whose pathologies are compatible.

In any case the builder declines every responsibility in case of use and working of the product in the following conditions, considered as **improper** use:

- displacement executed by non qualified staff;
- transport of those patients whose pathologies have been defined non compatible by a medical staff;
- non respect of the operators towards specifications in force about safety;
- transport of patients whose weight exceeds 170 Kgs;
- use of devices which may interfere with the stretcher's working;
- modifications or violations which has not been authorized by ultraMEDIC GmbH & Co. KG.
- non respect of what has been specified in this manual.

The scoop stretcher is essentially compound of:

- a main structure;
- an head zone profile;
- a bust zone profile;
- a feet zone profile.

Principal devices with relatives features are:

Locking of feet zone profile's extraction.

This device permits to adjust and lock the total length of the stretcher, in four positions, according to patient's height. The device is also utilized to get the stretcher folded, in stocking position.

Head side hook and feet side hook.

These devices permit the operators to open the stretcher in the required configuration (in two parts or in „A“ position) to place it under the patient without removing it, avoiding traumas which can come from his displacement.

Safety belts

These devices are utilized to limit voluntary or involuntary movements of the patient during the transport. The equipment includes three safety belts with quick-release, for a proper and safe immobilization of the patient.

2.3. TECHNICAL DATA

ultraSCOOP		
Height	mm	90
Width	mm	430
Minimum length	mm	1670
Maximum length	mm	2010
Folded length	mm	1200
Weight	kg	10
Maximum load	kg	170

These data are not binding for the builder, so they might be modified without notice.

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

2.2. TECHNICAL DESCRIPTION

The scoop stretcher, thanks to its particular features, has an important function of aid: in fact it permits an easy placement of the lifter under the patient without moving him and an easy immobilization during all transport. It is provided with devices, ergonomical and adjustable, which permit the hand-transport of patients always keeping a correct position, even in the most difficult conditions, as for example during the transport along staircases, tortuous and narrow ways or disconnected grounds. The correct position also consents to reduce operators' effort and so it brings to a situation of great security during the transport of the patient.

The structure made of anodized aluminium and the extremely strong construction guarantee an high grade of reliability.

1.5. RETURN OF MATERIALS

Before sending items of which it's required the renewal or repair under warranty, it is necessary a written authorization from Assistance Office of ultraMEDIC GmbH & Co. KG.

Faulty components must be properly packaged to avoid damaging during transport, sent ex works to our seat and provided with :

- Data drawn from the badge;
- Detailed description of fault and way in which it is happened.

The components for which warranty is valid are shipped ex works, the renewed ones remain a ultraMEDIC GmbH & Co. KG property.

1.6. EXCLUSIONS


Materials and components liable to ordinary wear and those whose lasting is not possible to determine before are not included in warranty.

Attention: Any form of warranty will decay if the badge is removed from the stretcher.

2. STRETCHER'S FEATURES AND TECHNICAL DATA

2.1. STRETCHER'S IDENTIFICATION

The stretcher is identified by an autoadhesive badge fixed on the main structure, which reports indelibly data about „CE“-marking.

The badge must never be removed. In case of damaging you must require a duplicate. ultraMEDIC GmbH & Co. KG doesn't recognize the stretcher if devoid of badge and also it let decay the warranty.

Returned goods / Replacements

Please check if all delivered parts are complete and functional.

ultraMEDIC GmbH & Co. KG is to be immediately informed of any damages or complaints.

Returned goods will only be accepted after approval with ultraMEDIC.

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
Use and maintenance instructions manual ultraSCOOP STRETCHER

1. WARRANTY

Warranty terms are stated in the following points, if not differently specified when executing the order.

1.2. THE OBJECT OF WARRANTY

The stretcher has been projected and built for a lasting utilization without particular problems; however in case of anomalies during the specified period covered by warranty, ultraMEDIC GmbH & Co. KG binds itself to repair or renew freely those parts which have presented break or premature wear caused by flaws in materials, fault during working, or faulting assembling. Warranty is not recognized for those parts whose break or wear have been caused by:

- non-respect of instructions contained in this use and maintenance manual;
- improper or undone maintenance;
- use of improper devices for ordinary and extraordinary maintenance;
- violations and modifications executed without ultraMEDIC GmbH & Co. KG authorization;
- utilization of non-original spare parts.

For commercial materials bought from external furnishers, ultraMEDIC GmbH & Co. KG recognizes to the buyer the same warranty as for ultraMEDIC GmbH & Co. KG suppliers.

1.3. DURATION OF WARRANTY

Warranty has a duration of 12 months and it starts 30 days after delivery day; the term is one and can't be delayed because of renewals or repairs executed during this period.

1.4. WAYS OF APPLICATION

To find the causes of anomalies and so to apply warranty the faulty components must be sent to ultraMEDIC GmbH & Co. KG. Repairs and renewals during warranty will be on ultraMEDIC GmbH & Co. KG opinion executed in its workshop or in assistance centers authorized by ultraMEDIC GmbH & Co. KG itself by thirds parties or in the same place. For operations executed in the same place the client will have to pay energy sources, extraordinary equipment and auxiliary staff.