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ONBO Electronic (Shenzhen) Co., Ltd.

No. 138, Huasheng Road, Langkou Community, Dalang Street, Longhua District, Shenzhen, China

Mircolife UAB, P. Lukšio g. 32, 08222 Vilnius, Lithuania

Please read these instructions carefully and completely before use and observe the care instructions.

Safety instructions

- The thermometer is suitable exclusively for measuring the body temperature!
- The minimum measurement time until the signaling tone sounds must be maintained without exception!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Therefore, never leave the thermometer within the reach of children.
- Protect the thermometer from shocks and impacts!
- Avoid ambient temperatures above 60°C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section “cleaning and disinfecting” to clean the thermometer.
- The thermometer must be intact when immersed in liquid disinfectant.
- Warning: The thermometer contains a long-life battery. For disposal, please comply with local regulations regarding hazardous waste.
- The manufacturer recommends having the accuracy checked by an authorized laboratory every 2 years.
- Do not store the unit in areas exposed to direct sunlight, high temperatures, high humidity, or dust. Performance may be negatively impacted.
- Using the unit in the immediate vicinity of mobile phones, microwave appliances, or other devices with strong electromagnetic fields may result in impaired functioning.
- Do not use the device if you think it is damaged or if anything appears unusual.

Turning on the thermometer

To turn on the thermometer, push the button over the display window; a short signaling tone means “thermometer ON”. A display test is then performed. At an ambient temperature of less than 32°C, an “l” and a flashing “°C” will appear on the display. The thermometer is now ready for use.

Usage

During a measurement, the current temperature is displayed continuously. The “°C” symbol flashes during this period. If the signaling tone sounds and the “°C” is no longer flashing, this means that the measured increase in temperature is now less than 0.1°C within 16 seconds and that the thermometer is ready for reading. To prolong the life of the battery, turn the thermometer off after use by pressing the operating button briefly. Otherwise, it will automatically turn itself off after approximately 9 minutes.

Storing measured values

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small “M” for memory will appear at the right of the display. After 3 more seconds, this value will disappear and the thermometer will switch to the normal measurement mode. The stored value will then be erased.

Cleaning and disinfecting

Name: Isopropyl alcohol 70%

Use a cotton swab moistened with alcohol (70% isopropyl) to clean and disinfect the thermometer casing and the measuring probe. Never use abrasive detergents, thinners, or benzene for cleaning and **never immerse the instrument in water or other cleaning liquids.**

Replacing the battery

When the “q” symbol (upside-down triangle) appears in the lower right of the display, the battery is empty and needs replacing. To replace the battery, please remove the compartment lid of the thermometer. Replace the empty battery by positioning the plus at the top. Make sure a battery of the same type is available. Batteries can be purchased at any electrical goods shop. For battery disposal, please refer to the safety instructions..

Guarantee

The thermometer is guaranteed for two years from the date of purchase. The guarantee does not apply to damage caused by improper handling. The batteries and the packaging are also excluded from the guarantee. All other damage claims are excluded. Any guarantee claim must be submitted with the purchase receipt.

Technical data

Art. no.:	MT16F1/5810090
Type:	Maximum thermometer
Measurement range:	32°C to 43.9°C Temperature < 32.0°C: Display “l” for low (too low) Temperature > 43.9°C: Display “H” for high (too high)
Measurement accuracy:	±0.1°C between 32.0°C and 43.9°C
Self test:	Automatic internal test of test value 37°C; if there is a deviation of > 0.1°C, “ERR” (error) is displayed.
Display:	LCD display with 3 digits: lowest display unit 0.1°C
Signaling tone:	To test operational readiness and if the temperature rises by less than 0.1°C within 16 seconds
Memory:	To automatically store the last measured value
Operating conditions:	+16°C to +40°C, 15-95% maximum relative air humidity
Storage temperature:	-25°C to +60°C, 15-95% maximum relative air humidity
Battery:	1.5/1.55 V; lifetime: Approx. 4,500 measurements
Reference to standards:	IEC 60601-1; IEC 60601-1-2, IEC 60601-1-11 Electromagnetic compatibility: Device fulfills the stipulations of the standard IEC 60601-1-2, EN 12470-3, clinical thermometers, ASTM standard

The device has an expected lifetime of 5 years and complies with the stipulations of the EU Directive 93/42/EEC.

Westmed distributed by Fleischhacker GmbH & Co. KG

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