



WESTMED

DATA SHEET WESTMED® Rappaport Stethoscope



Order number

5810026

Description

Our WESTMED® stethoscope type Rappaport with high-quality double tube system guarantees optimum sound conduction. To switch the funnel, the double-head chest-piece is rotated by 180°. In general, this chest piece is used to detect low frequency heart sounds, diastolic and third heart sounds as well as first and second heart sounds. The high-pitched heart sounds emanating from the hip can also be detected. The stethoscope is supplied with three sound funnels in small, medium and large sizes. This gives you flexibility in use depending on the age of the patient. Heart sounds can be listened to from infants to adults.

- Durable tubes with metal clip

| Technical information | |
|----------------------------------|---|
| Reference number | HS-30C1 |
| Color | Black |
| Total length | 76-77 cm |
| Ambient temperature and humidity | +10 °C to +40 °C; 10% to 90% |
| Storage temperature and humidity | -20 °C to +70 °C; 10% to 90% |
| Durability | More than 10,000 measurements or 5 years from production date |
| Packaging unit | 1 pieces per pack 50 packs per carton 10 cartons per pallet |



WESTMED

DATA SHEET WESTMED® Rappaport Stethoscope

| Receiver | |
|---------------------|-------------------|
| Material | Anodized aluminum |
| Weight | 37.5 g |
| Diameter large side | 47 mm |
| Diameter small side | 30 mm |

| Multi-frequency diaphragm | |
|---------------------------|---------|
| Material | PVC |
| Diameter large side | 45.5 mm |
| Diameter small side | 31.5 mm |

| Cold protection rings | |
|-----------------------|-------------------------|
| Material | Thermo-propylene rubber |
| Diameter large side | 47.3 mm |
| Diameter small side | 37.5 mm |

| Remaining components | |
|----------------------|--------|
| Length of fork | 60 cm |
| Material fork | Copper |
| Material feather | Copper |
| Material tube | PVC |

| General information | |
|---------------------|--|
| Safety notes | Do not immerse the stethoscope in any liquid. Do not use the device if it is damaged. The device cannot be sterilized. |









| Regulatory information | |
|------------------------|---|
| Classification | Class I medical device According to Annex VIII, Regulation (EU) 2017/745 |

| Packaging symbols | |
|-------------------|----------------------------------|
| MD | Medical device |
| REF | Reference number/article number |
| EC REP | Authorized representative Europe |



WESTMED

DATA SHEET WESTMED® Rappaport Stethoscope

| | |
|--|--|
|  | Distributor |
|  | Importer |
|  | Manufacturer |
|  | Date of production |
|  | Observe the operating instructions |
|  | Protect against moisture |
|  | Protect against heat |
|  | The CE marking is short for „Conformité Européenne“, which means „European Conformity“. It symbolizes the conformity of the product with the applicable requirements that the European Community places on the manufacturer. |



HONSUN (NANTONG) Co.,
Ltd No. 8 Tongxing Road Economic & Technological
Development Area, 226009 Nantong City, Jiangsu Province, P. R. China



Shanghai International Holdin Corp. GmbH (Europe)
Eiffestraße 80, Hamburg, Germany