

Model FED 115 | Drying and heating chambers with forced convection and enhanced timer functions

BINDER heating chambers of the FED Avantgarde.Line series provide almost unlimited capacity and are particularly adaptable. With advanced timing functions and a controllable fan, temperature and convection conditions are easily controlled.

BENEFITS

- · Best energy efficiency on the market
- · Excellent temporal and spatial temperature accuracy
- USB connection for recording data
- Temperature range up to 300 °C



Model 115

MAIN FEATURES

- Temperature range: ambient temperature plus 10 °C to 300 °C
- Up to 30% lower energy consumption compared to conventional units on the market
- High temperature accuracy thanks to APT.line[™] technology
- · Adjustable fan speed
- · Controller with LCD display and enhanced time functions
- Electromechanical control of the exhaust air flap
- · 2 chrome-plated racks

- Stackable
- Class 2 integrated independent adjustable temperature safety device
 (DIN 12880) with visual alarm
- · Ergonomic handle design
- Ethernet interface
- · USB port for recording data

ORDERING INFORMATION

| Interior volume [L] | Voltage | Option model | Version | ArtNo. |
|---------------------|----------------------|--------------|---------------|-----------|
| 116 | 230 V 1~ ph 50/60 Hz | Standard | FED115-230V | 9010-0293 |
| 110 | 120 V 1~ ph 60 Hz | Standard | FED115UL-120V | 9010-0294 |




TECHNICAL DATA

| Description | FED115-230V ¹ | FED115UL-120V ¹ |
|---|--------------------------|----------------------------|
| Article Number | 9010-0293 | 9010-0294 |
| Performance Data Temperature | | |
| Temperature range 10 °C above ambient temperature to [°C] | 300 | 300 |
| Temperature variation at 150 °C [± K] | 1.2 | 1.2 |
| Temperature fluctuation at 150 °C [± K] | 0.3 | 0.3 |
| Heating-up time to 150 °C [min] | 18 | 18 |
| Recovery time after 30 seconds door open at 150 °C [min] | 4 | 4 |
| Air change data | | |
| Air change (approx.) at 100 °C [x/h] | 32 | 32 |
| Electrical data | | |
| Rated Voltage [V] | 230 | 120 |
| Power frequency [Hz] | 50/60 | 60 |
| Nominal power [kW] | 1.3 | 1.4 |
| Unit fuse [A] | 6.3 | 12.5 |
| Phase (Nominal voltage) | 1~ | 1~ |
| Measures | | |
| Interior volume [L] | 116 | 116 |
| Net weight of the unit (empty) [kg] | 57 | 57 |
| Permitted load [kg] | 150 | 150 |
| Load per rack [kg] | 30 | 30 |
| Wall clearance back [mm] | 160 | 160 |
| Wall clearance sidewise [mm] | 100 | 100 |
| Internal Dimensions | | |
| Width [mm] | 550 | 550 |
| Height [mm] | 550 | 550 |
| Depth [mm] | 385 | 385 |
| Doors | | |
| Unit doors | 1 | 1 |
| Housing dimensions not incl. fittings and connections | | |
| Width net [mm] | 710 | 710 |
| Height net [mm] | 705 | 705 |
| Depth net [mm] | 605 | 605 |
| Environment-specific data | | |
| Energy consumption at 150 °C [Wh/h] | 340 | 340 |
| Sound-pressure level [dB(A)] | 43 | 43 |
| Fixtures | | |
| Number of shelves (std./max.) | 2/5 | 2/5 |

1 All technical data is specified for unloaded units with standard equipment at an ambient temperature of +22 °C ±3 °C and a power supply voltage fluctuation of ±10 %. The temperature data is determined to BINDER factory standard following DIN 12880, observing the recommended wall clearances of 10 % of the height, width, and depth of the inner chamber. Technical data refers to 100 % fan speed. All indications are average values, typical for units produced in series. We reserve the right to change technical specifications at any time.



DIMENSIONS incl. fittings and connections [mm]



OPTIONS

| Designation | Description | * | ArtNo. |
|---|--|--------|-----------|
| | left | | |
| | 10 mm | 01 | 8012-1278 |
| | 30 mm | 01 | 8012-1655 |
| | 50 mm | 01 | 8012-1657 |
| | 100 mm | 01 | 8012-1669 |
| | right | | |
| | 10 mm | 01 | 8012-1271 |
| Access port with silicone plug | 30 mm | 01 | 8012-1654 |
| | 50 mm | 01 | 8012-1656 |
| | 100 mm | 01 | 8012-1670 |
| | top | | |
| | 10 mm | 01 | 8012-1264 |
| | 30 mm | 01 | 8012-1365 |
| | 50 mm | 01 | 8012-1381 |
| | 100 mm | 01, 10 | 8012-1672 |
| Alarm function for overheating | Switchable acoustic alarm, with adjustable limit value on the independent temperature safety device | - | 8012-1637 |
| Analog output 4-20 mA | for temperature values (output not adjustable) | 02 | 8012-1625 |
| Calibration certificate, expanded | for temperature; for extending the measurement in center of chamber to include another test temperature | - | 8012-1118 |
| | for temperature, measurement in center of chamber at specified temperature | - | 8012-1137 |
| Calibration certificate, | temperature measurement incl. certificate and 27 measuring points at specified temperature | - | 8012-1596 |
| temperature | temperature measurement incl. certificate, 15-18 measuring points at specified temperature | - | 8012-1575 |
| | temperature measurement incl. certificate, 9 measuring points at specified temperature | - | 8012-1554 |
| Class 3.1 independent temperature safety device | with visual alarm (DIN 12880) | 04, 12 | 8012-1644 |
| Clock | battery-backed | - | 8012-1642 |
| Door gasket | made of FKM, silicone-free | 03 | 8012-1651 |
| Door lock | lockable door handle | - | 8012-1666 |
| Gas-tight option model | Incl. additional measures for the greatest possible impermeability to gas; not in conjunction with access port or interior lighting option | 09 | 8012-1686 |
| Heating output, | For faster heating up time and improved recovery time; max. 300 °C; not in conjunction with load management option | | |
| reinforced | 120/230 V option model | - | 8012-1699 |
| HEPA air filter | On fresh air supply for unit; filter classification H14 (acc. to. EN 1822-1:2009, > 99.995% at 0.3 μm) | - | 8012-1696 |
| Inert gas connection | With gas inlet and outlet, Ø 10 mm, incl. additional measures for the greatest possible impermeability to gas; not in conjunction with access port or interior lighting option | - | 8012-1692 |
| nner chamber, reinforced | max. total load 170 kg | - | 8012-1714 |
| | | | |

* Notes > See last page

BINDER

| Designation | Description | * | ArtNo. |
|--------------------------------------|--|---|-----------|
| | with 15 W light bulb | · | |
| Interior lighting | 120 V option model | - | 8012-1633 |
| | 230/240 V option model | - | 8012-1629 |
| Measurement of air exchange rate | in accordance with ASTM D5374, definition and protocol according to ambient temperature | - | 8012-1197 |
| Pt 100 temperature sensor | additional flexible Pt 100, interior, for displaying the temperature on the unit display | - | 8012-1621 |
| Viewing window and interior lighting | viewing window (triple insulating glass) in door, 370 x 370 mm, and 15 W interior | | |
| | 120 V option model | - | 8012-1707 |
| | 230 V option model | - | 8012-1705 |

* Notes > See last page

ACCESSORIES

| Designation | Description | * | ArtNo. |
|--------------------------|---|----|----------|
| APT-COM™ 4 BASIC- | for simple logging and documentation requirements with up to 5 networked units. | | |
| Edition | version 4, BASIC edition | 19 | 9053-003 |
| APT-COM™ 4 GLP- | for working under GLP-compliant conditions. Measured values are documented in a tamper-proof way in line with the requirements of FDA Regulation 21 CFR 11. | | |
| Edition | version 4, GLP edition | 19 | 9053-004 |
| APT-COM™ 4 | convenient unit and user management built on the BASIC edition. Suitable for networking up to 100 units. | | |
| PROFESSIONAL- Edition | version 4, PROFESSIONAL edition | 19 | 9053-004 |
| Data Logger Kit | T 350: For continuous temperature logging from 0 °C to 350 °C; kit includes 1 data logger, Pt 100 sensor with 2 m extension cable and 1 magnetic fixture for mounting to the BINDER unit | 19 | 8012-071 |
| Data Logger Software | LOG ANALYZE software kit, configuration and evaluation software for all BINDER Data Logger Kits (incl. USB data cable) | 19 | 8012-082 |
| | with lid made from 18/8 stainless steel, stackable, dimensions: 230 x 133 x 45 mm | - | 4022-012 |
| Instrument tray | with lid made from 18/8 stainless steel, stackable, dimensions: 315 x 215 x 60 mm | - | 4022-012 |
| pH-neutral detergent | concentrated, for gentle remove of residual contaminants; 1 kg | - | 1002-001 |
| Qualification documents | IQ/OQ documents – supporting documents for validation performed by customers, consisting of: IQ/OQ checklists incl. calibration guide and comprehensive unit documentation; parameters: temperature, CO ₂ , O ₂ , pressure, depending on unit | | |
| | Digital in PDF format | - | 7057-000 |
| | Hard copy inside folder | - | 7007-000 |
| Deek | chrome plated | - | 6004-016 |
| Rack | stainless steel | - | 6004-015 |
| Rack, heavy load | Stainless steel, max. load per rack 70 kg | - | 6004-019 |
| Rubber pads | set anti-slip feet | - | 8012-188 |
| Shelf, perforated | Stainless steel | - | 6004-018 |
| Shell, periorated | | | |

* Notes > See last page

SERVICES

| Description | * | ArtNo. |
|--|---|--|
| | | |
| Unpacking and setting up of unit, connecting to existing connections, and checking function | 13, 18 | DL10-0100 |
| Unit function instructions for operation and programming of the controller | 18 | DL10-0500 |
| | | |
| Functional testing of all electrical and mechanical components, short calibration, documentation in the maintenance schedule | 14, 18 | DL20-0200 |
| | | |
| including certificate (in accordance with ASTM D5374) | 14, 16, 17, 18 | DL33-0000 |
| | Unpacking and setting up of unit, connecting to existing connections, and checking function Unit function instructions for operation and programming of the controller Functional testing of all electrical and mechanical components, short calibration, documentation in the maintenance schedule | Description 13, 18 Unpacking and setting up of unit, connecting to existing connections, and checking function 13, 18 Unit function instructions for operation and programming of the controller 18 Functional testing of all electrical and mechanical components, short calibration, documentation in the maintenance schedule 14, 18 including certificate (in accordance with ASTM D5374) 14, 16, |

BINDER

| Designation | Description | * | ArtNo. |
|---|--|-------------------|-----------|
| - | Expansion – including certificate, each additional measuring point in center of chamber at specified test temperature | 14, 16, 17, 18 | DL30-0102 |
| Temperature calibration | Temperature calibration with 1 measuring point in center of chamber at a specified test temperature, including certificate | 14, 16, 17, 18 | DL30-0101 |
| Temperature measurement, 18 measuring points | Temperature measurement with 18 measuring points at a specified test temperature, including certificate | 14, 16, 17, 18 | DL30-0118 |
| Temperature measurement, 27 measuring points | Temperature measurement with 27 measuring points at a specified test temperature, including certificate | 14, 16, 17, 18 | DL30-0127 |
| Temperature measurement, 9 measuring points | Temperature measurement with 9 measuring points at a specified test temperature, including certificate | 14, 16, 17, 18 | DL30-0109 |
| Validation services | | | |
| Execution of IQ/OQ | Execution of IQ/OQ in accordance with qualification folder | 15, 18 | DL40-0100 |
| Execution of IQ/OQ/PQ | Execution of IQ/OQ/PQ in accordance with qualification folder | 15, 18 | DL44-0500 |
| Warranty service | | | |
| Warranty extension from 2 to 3 years | The warranty is extended from 2 to 3 years from the delivery date, wear parts are excluded | - | DL00-4041 |
| Warranty extension from 2 to 5 years | The warranty is extended from 2 to 5 years from the delivery date, wear parts are excluded | - | DL00-4042 |
| | | | |

* Notes > See last page



NOTES

- Condensation may occur in the area around the access port. Access ports may be placed in custom locations for an additional charge. 01
- 02 UL mark is not granted when this option is used.
- Heat resistant only to max. 200 °C. 03
- Only available on units rated for 230 V. 04
- 07 The additional heat input may affect the temperature behavior. Not in conjunction with the optional access port, door with window and interior lighting.
- 09 10 Not available on 23-liter units.
- Not available on 23- or 53-liter units. Only available on units rated for 230 V or 400 V. 11 12
- 13 Installation and connections take place at unit location; transport within the company only upon consultation.
- 14 We recommend a BINDER service contract to cover unit inspections, calibrations and validations.
- 15 OQ according to Yellow Paper = completed factory validation documentation of all OQ checklists.
- 16 Sensor calibration is performed in an accredited calibration laboratory. Calibration is performed according to the BINDER factory standard. 17
- Quoted prices do not include travel costs. Please refer to the chapter on BINDER Service for travel costs for your region. Quoted prices for services performed in Switzerland do not 18 include a country-specific added fee (available on request).
- For additional accessories, refer to the Process documentation chapter. 19

BINDER GmbH Tuttlingen, Germany TEL +49 7462 2005 0 FAX +49 7462 2005 100 info@binder-world.com www.binder-world.com

BINDER Asia Pacific (Hong Kong) Ltd. Kowloon, Hong Kong, P.R. China TEL +852 39070500 FAX +852 39070507 asia@binder-world.com www.binder-world.com

BINDER Environmental Testing Equipment (Shanghai) Co., Ltd. Shanghai, P.R. China TEL +86 21 685 808 25 FAX +86 21 685 808 29 china@binder-world.com www.binder-world.com

Representative Office for CIS Moscow. Russia TEL +7 495 988 15 16 FAX +7 495 988 15 17 russia@binder-world.com www.binder-world.com

BINDER Inc. Bohemia, NY, USA TEL +1 631 224 4340 FAX +1 631 224 4354 usa@binder-world.com www.binder-world.us