



# Hunting is over!

Introducing our Space Savers




reddot winner




# STAVANGER Ampoule kit

## In short: Organised

Our new STAVANGER ampoule systems inspire with their clear arrangement and organization of medication and all other materials required. The unique retaining clips developed by us protect the sensitive glass ampoules and guarantee quick removal in an emergency case. The clips can be arranged both lengthwise and crosswise to make optimal use of the space available. Never before has it been so easy to store many ampoules. With the support of labeling clips, refilling is child's play. The ampoule arsenals are available in different sizes with different fittings. Your advantage: All individual parts can be reordered.



LoopFlex rubber bands for fixing larger utensils



Various foil pockets for the storage of loose materials

SingleFlex rubber bands in different sizes for the fixation of other kinds of ampoules

Transparent pockets in different sizes for small items

**STAVANGER Ampoule Kit XL**  
up to 300 Ampoules  
L 25 x W 36 x H 12 cm

**STAVANGER Ampoule Kit MM**  
up to 150 Ampoules  
L 18 x W 24 x H 12 cm

Quick overview using labelling clips with ingredient stickers

Dip and wipe disinfection of all components possible

Optimum use of space thanks to the adjustable retaining clips which can be used lengthwise and crosswise

Optional middle wall for additional space

Removable system plate for hygienic reception in the magnetic system

Robust distance clips to prevent glass breakage

Complete, thermo-coated outer shell





FLEISCHHACKER

**Fleischhacker GmbH & Co. KG**

An der Silberkuhle 18  
58239 Schwerte

T 02304 931-0  
F 02304 931-199

stavanger@fleischhacker.biz  
www.fleischhacker.biz



**STAVANGER V4000**

System bag for variable use, for example for drugs that require intubation equipment or child emergency  
W 52 x H 37 x D 12 cm

**Interested?**

Book your online presentation now.

<https://www.fleischhacker.biz/de/online-presentations>

